

ENHANCE PATIENT SAFETY DURING UPPER GI PROCEDURES



Microstream™ Advance bite block filter line

Every patient is unique. And some patients may be more sensitive to sedation during upper gastrointestinal procedures. That's why we offer Microstream™ Advance bite block filter line.



The Microstream™ Advance bite block filter line is designed to:

- Capture etCO₂ samples from monitors enabled with Microstream™ capnography.
- Deliver supplemental oxygen.
- Help protect endoscopes.
- Maintain stability without compromising endoscope maneuverability.

Features:

- Provides continuous CO₂ monitoring to alert you to changes in adequacy of ventilation.
- Integrated O₂ tubing for supplemental oxygen offers high flow delivery (>10 liters per minute).†
- O₂ tubing offered in standard length (6.5 ft/200 cm) and long-length (13 ft/400 cm).
- Accepts endoscopes up to 60 French.

† If O₂ is delivered at flow rates greater than 10 lpm, the patient's exhaled CO₂ may be diluted, producing attenuated waveforms and etCO₂ values.

© 2021 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 21-weu-guardian-ad-4986092

Medtronic
Further, Together